



# Eye tracking


Woldorff lab meeting  
October 26, 2010

# Reasons to use eye-tracking

- Ensure subject maintains fixation
  - Mean eye-position, std dev of eye-position
  - Time-locked average eye-trace to left-cue against right-cue
- Follow saccades to regions of interest (ROIs)
  - Calculate ROI dwell time
- Generate “heat maps”
- Generate “overlay traces”

# Planning covert shifts of attn.

EOG  
Data



# Traces as a function of task



Free examination.

1



2

Estimate material circumstances of the family



3

Give the ages of the people.



4

Surmise what the family had been doing before the arrival of the unexpected visitor.



5

Remember the clothes worn by the people.



6

Remember positions of people and objects in the room.



7

Estimate how long the visitor had been away from the family.

3 min. recordings of the same subject

# Visual processing and disorders

Face scanning  
(autism group)



Face scanning  
(control group)





# Mapping salience




<http://www.etre.com/usability/eyetracking/showme/>

# Biological salience

women




men



# BIAC Eyetracking Equipment


- Camera w/IR LED illumination mounted to head-coil (30-60Hz scanning, adj.)
- NTSC signal (interlaced) sent to TV-capture board
- Viewpoint interface uses video signal to segment pupil and track its location
- After calibration, VP records location data to file.





<http://fourier.biac.duke.edu/wiki/doku.php/biac:experimentalcontrol:eyetrackviewpoint>


# Step 1: get good pupil contrast



Make sure that the eye is well-centered in the visual field. This needs to be done before the subject is slid into the scanner, and the MR tech will help you with this part. Then, use the console to adjust the LED intensity if necessary.

Toggle “show threshold dots” and adjust “Brightness” and “Contrast” sliders to segment the pupil. Adjust bounding box to constrain search area.

# Step 2: Calibrate well



Instruct subject on “following the blinking cues”. Warn subject that calibration is about to start. Do twice in a row so that subject is familiar with procedure. Ensure that the calibration matrix that comes out is regular-looking and not folded back upon itself. (example depicts a decent calibration). Try adjusting parameters of calibration routine.

REPEAT AS MANY TIMES AS NECESSARY UNTIL YOU ARE SATISFIED!!

## Step 3:

- Follow instructions on the BIAC wiki to get Presentation to communicate with Viewpoint so that stimulus events can be recorded in the eye-tracking datastream.
- Double-check timing, and be sure that you can keep both the Presentation and the Viewpoint logfiles in sync!
- Check logfile for markers after running
- Keep an eye on traces
- (this is harder than you might think)

# Problems:

- Head movement throws off calibration
- Hard to double-check afterward, as the possibly-faulty coordinates are the only things recorded (and not the video signal).
- Not so great software
  - Crashy? Please record error message and give to me.
  - from BIAC talk, JV's calibration routines superior to VP's
  - (for \$7K, frankly, I'm not impressed.)
- Positioning: only MR tech, fixed at too-great distance
- Not done with so much frequency.



# Drift correction

- Recalibrate every run?
- High-pass filter correction on y-coordinates?
- Adjust the calibration mtx at every fixation?

“True” data:



Corrupted data:



Trend-removal:



# Software we have:

- Presentation + VP (codes to VP over serial)
  - Somewhat difficult, what we do now.
- CIGAL + VP (eye locs over serial to CIGAL)
  - Can do calib, replay traces on actual pdigm, display horz+vert displacement w/physiological data and MR trigs
- Presentation + VP-Pres Addon + VP
  - Easy integration w/ presentation, not installed yet
- Chris Petty's script for ROI analysis, bxh\_events
- Ken's scripts for quick MATLAB analysis (unfinished)


# CIGAL Replay



# CIGAL Replay




# Define an ROI (Chris P.)





# Ruth's SaccAtt Experiment




# Viewpoint logfile (txt or xls)


	A	B	C	D	E	F	G	H	I	J
1	3	Product Version: PC60								
2	3	Executable File Version: 2.6.1								
3	3	Program Build Date: Aug 13 2002, 20:41:02								
4	3	Serial Number: 200-326-2-2								
5	3	Customer Name: Jim Voyvodic, Duke University, NC								
6	3	-----								
7	3	TimeValue	2008	10	23	14	4	56	UTC	
8	3	TimeStam	Thursday, October 23, 2008, 10:04:56 AM							
9	3	DataForma	57							
10	3	-----								
11	3	TotalTime	DeltaTime	X_Gaze	Y_Gaze	Region	PupilWidth	PupilAspe	Count	Torsion
12	2	0 +								
13	10	0.0147	0	0.4735	0.484	0	0.1562	0.9615	0	-998
14	10	0.048	33.263	0.4709	0.3515	-1	0.1625	1	1	-998
15	10	0.0815	33.492	0.4709	0.3515	-1	0.1562	0.9259	2	-998
16	10	0.1147	33.22	0.4529	0.3216	-1	0.1594	0.9273	3	-998
17	10	0.148	33.351	0.4691	0.2876	-1	0.1562	0.9259	4	-998
18	10	0.1815	33.418	0.4358	0.324	-1	0.1625	0.9811	5	-998
19	10	0.2143	32.886	0.4538	0.3538	-1	0.1594	0.9444	6	-998
20	10	0.248	33.679	0.4717	0.384	-1	0.1562	0.9091	7	-998

# Lost pupil and markers


42302	10	1398.509	33.702	1.8685	-0.6242	-1	0.0312	0.8	41117	-998
42303	99	1398.541	3	0						
42304	10	1398.575	66.441	2.2421	-0.5407	-1	0.0406	0.6154	41118	-998
42305	10	1398.608	33.212	0.4306	-0.0864	-1	0.0969	0.6129	41119	-998
42306	10	1398.642	33.446	1.7151	-0.5053	-1	0.0094	0.3333	41120	-998
42307	10	1398.675	33.346	1.875	-0.254	-1	0.0281	0.1111	41121	-998
42308	99	1398.708	3	0						
42309	99	1398.741	3	0						
42310	99	1398.775	3	0						



42534	10	1406.205	33.518	1.0986	0.1939	-1	0.1969	0.8095	41322	-998
42535	12	1406.221	1							
42536	10	1406.238	33.125	1.1551	0.1905	-1	0.2125	0.75	41323	-998
42537	10	1406.272	33.617	0.7679	0.2526	-1	0.1469	0.9038	41324	-998




# est\_quality




Quality information for file: 35104\_1to4.xls

Number of points:	58007
Points with 'no data':	1514
Glitches:	2
Events:	567
Gaps > two samples:	90
Percent <0 and >1:	
x:	3.48, 5.34
y:	0.87, 39.51
Percent <-0.5 and >1.5:	
x:	0.37, 1.30
y:	0.10, 12.53




# est\_quality




Quality information for file: 35706\_1to4.xls

Number of points:	57780	
Points with 'no data':	238	
Glitches:	1	
Events:	563	
Sample intervals > 2*period:	1	
Percent <0 and >1:		
x:	15.98,	10.84
y:	46.49,	7.90
Percent <-0.5 and >1.5:		
x:	3.16,	0.06
y:	6.68,	5.26







# plot\_trace




# plot\_traces




# plot\_x\_histograms



# plot\_joint\_histograms



# Make\_pst





# Make\_pst

